

1CCD PORTABLE USB CAMERA



Video to USB



USB CAMERA Technical Data

Video system: PAL/NTSC
Image sensor: 1/3"CCD
Horizontal resolution: 480 lines
Min illumination: 0.4Lux
Aspect ratio: 4:3
SNR: ≥56dB
White balance: Auto/Manual
Gain: AUTO/OFF
Shutter speed: 1/60S-1/10,000S
Output : Video ×1

CAMERA HEAD

Dimension	Φ30×80mm
Weight	230g
Lens connection	C-mount
Cable length	3m
Protection	IPX8

CONTROL UNIT

Dimension	115mm×180mm×48mm
Weight	530g
Protection class	IP20
Line voltage	220V/50Hz
Power consumption	220V/50Hz
Power tested acc. To	UL2601-1, en60601-1, IEC601-1-2



Clarity In Detail-Brilliance In View

SINCE 1998

HD PORTABLE USB CAMERA



Video to USB

USB HD CAMERA Technical Data

Video system: PAL/NTSC
Image sensor: 1/3"CMOS
Horizontal resolution: 1000 lines
Min illumination: 0.4Lux
Aspect ratio: 16:9
SNR: ≥56dB
White balance: Auto/Manual
Gain: AUTO/OFF
Shutter speed: 1/60S-1/10,000S
Output : HDMI ×1 VIDEO X1

CAMERA HEAD

Dimension	Φ30×80mm
Weight	230g
Lens connection	C-mount
Cable length	3m
Protection	IPX8

CONTROL UNIT


Dimension	115mm×180mm×48mm
Weight	530g
Protection class	IP20
Line voltage	220V/50Hz
Power consumption	220V/50Hz
Power tested acc. To	UI2601-1, en60601-1, IEC601-1-2



Clarity In Detail-Brilliance In View

SINCE 1998

1CCD & HD PORTABLE CAMERA



1CCD Technical Data

Video system : PAL/NTSC
Image sensor : 1/3"CCD
Horizontal resolution : 700 lines
Min illumination : 0.01Lux
Aspect ratio : 4:3
SNR : ≥48dB
White balance : Auto/Manual
Gain : AUTO/OFF
Shutter speed : 1/60S-1/10,000S
Video outputs : VIDEO×1
Temperature : -20℃ ~ 60℃
Humidity : ≤80%

HD Technical Data

Video system : PAL/NTSC
Image sensor : CMOS
Horizontal resolution : 1920X1080
Min illumination : 0.017Lux
Aspect ratio : 16:9
SNR : ≥48dB
White balance : Auto/Manual
Gain : AUTO/OFF
Shutter speed : 1/60S-1/10,000S
Video outputs : VIDEO×1 HD-SDI X1
Temperature : -20℃ ~ 60℃
Humidity : ≤80%



Clarity In Detail-Brilliance In View

SINCE 1998

7

Beijing Fanxing Guangdian Medical Treatment Equipment Co.,Ltd.

86-1358-177-6026

DELONDONG@HOTMAIL.COM



86-1358-177-6026



DELONDONG